

MIRS Measurement on Cassava Cell Walls and Flours

High-Throughput Phenotyping Protocols (HTPP), WP3

Montpellier, France, November 2020

Karima MEGHAR, Centre de Coopération Internationale en Recherche Agronomique pour le Développement (CIRAD), Montpellier, France



www.rtbfoods.cirad.fr

cirad

This report has been written in the framework of RTBfoods project. (Arial 12, lowercase)

To be cited as:

Karima MEGHAR, (2021). *MIRS Measurement on Cassava Cell Walls and Flours. High-Throughput Phenotyping Protocols (HTPP), WP3*. Montpellier, France: RTBfoods Laboratory Standard Operating Procedure, 12 p.

Ethics: The activities, which led to the production of this document, were assessed and approved by the CIRAD Ethics Committee (H2020 ethics self-assessment procedure). When relevant, samples were prepared according to good hygiene and manufacturing practices. When external participants were involved in an activity, they were priorly informed about the objective of the activity and explained that their participation was entirely voluntary, that they could stop the interview at any point and that their responses would be anonymous and securely stored by the research team for research purposes. Written consent (signature) was systematically sought from sensory panelists and from consumers participating in activities.

Acknowledgments: This work was supported by the RTBfoods project <https://rtbfoods.cirad.fr>, through a grant OPP1178942: Breeding RTB products for end user preferences (RTBfoods), to the French Agricultural Research Centre for International Development (CIRAD), Montpellier, France, by the Bill & Melinda Gates Foundation (BMGF).

Image cover page © LAJOUS P. for RTBfoods.

RTBfoods WP3: High-Throughput Phenotyping Protocols (HTPP)		
SOP: MIRS Measurement on Cassava Cell Walls and Flours		
Date: 04/11/2020		Release: 1
Written by Karima MEGHAR For information on this SOP please contact: Karima MEGHAR, karima.meghar@cirad.fr		
This document has been reviewed by:		
Fabrice DAVRIEUX (CIRAD) Karima MEGHAR (CIRAD) Emmanuel ALAMU (IITA) Fabrice DAVRIEUX (CIRAD)		06/11/2020 09/11/2020 13/11/2020 17/11/2020
Final validation by:		
Fabrice DAVRIEUX (CIRAD)		17/11/2020

CONTENTS

Table of contents

1	Scope and application.....	6
2	References	6
3	Principle.....	7
4	Procedure	7
4.1	MIRS measurements.....	7
4.2	Repeatability.....	8
5	Critical points or Note on the procedure.....	9
6	Appendices	10
6.1	Appendix 1: Procedure of extraction of cell walls from cassava flour	10

ABSTRACT

The purpose of this SOP is to use MIRS to predict the cooking behaviour of cassava from different varieties by using spectra of corresponding cassava cell walls or flour. For this purpose, eighty-nine samples of cell walls from 37 different varieties (2 or 3 roots per genotype) provided by CIAT (Colombia) were analysed by using Thermo Scientific™ Nicolet™ iS50R Research FTIR Spectrometers in the 4000-400 cm⁻¹ region.

For the repeatability test, 8 spectra of the same flour sample were acquired in order to determine the optimum number of MIRS measurements (replicates) required to be representative of the sample. The RMS obtained for 5 combinations of two replicates (2-3, 1-2, 3-5, 3-7 and 6-8) are 457, 1066, 710, 1352 and 614 μ abs respectively. The RMS values between 2 replicates, randomly selected, were lower than the RMS mean for all the replicates. These results indicate that two replicates of the sample result to a spectral dispersion (variability) similar to the one obtained with 8 repetitions.

Key words: MIRS, Cassava cell walls, Cooking behaviour.

1 SCOPE AND APPLICATION

One of the fingerprinting techniques widely used to characterize food ingredients and detect possible adulterants is Fourier transform infrared spectroscopy (FT-IR) in the middle range (400–4000 cm⁻¹, FTMIR). FT-MIR monitors the fundamental vibrational and rotational stretching modes of molecular bonds, which produce a chemical profile of the sample and provide a greater amount of chemical information than FT-NIR (Lohumi, Lee, Lee, & Cho, 2015). Thus, FT-MIR has been used to classify cereal flours, such as wheat, oats and buckwheat (Cocchi et al., 2004), to discriminate different types of wheat flours from rye and triticale flours (Sujka, Koczoń, Ceglińska, Reder, & Ciemniowska-Zytkiewicz, 2017) or to determine components of quinoa flour (González-Muñoz et al., 2016). Trans-fatty acids in cereal products were determined by (Kim, Himmelsbach, et al, 2007). Suchowilska, Kandler, Wiwart, and Krska (2012) could differentiate between species of wheat by applying Principle Components Analysis to MIR spectra. Tamaki and Mazza (2011) quantified different types of polysaccharides, ash, and extractives in straw, albeit with mixed results.

The purpose of this SOP is to use MIRS to predict the cooking behaviour of cassava from different varieties by using MIRS spectra of corresponding cassava cell walls or flour. This SOP describes the MIR measurement on Cassava cell wall and flour and how can we discriminate between good and bad cooking time by using MIR spectra data. In this protocol, cell walls analysis are taken as an example.

2 REFERENCES

- Cocchi, M.; Foca, G.; Lucisano, M.; Marchetti, A.; Pagani, M. A.; Tassi, L.; Ulrici, A. Classification of Cereal Flours by Chemometric Analysis of MIR Spectra. *J. Agric. Food Chem.* 2004, 52 (5), 1062–1067. <https://doi.org/10.1021/jf034441o>.
- González-Muñoz, A.; Montero, B.; Enrione, J.; Matiacevich, S. Rapid Prediction of Moisture Content of Quinoa (*Chenopodium Quinoa* Willd.) Flour by Fourier Transform Infrared (FTIR) Spectroscopy. *Journal of Cereal Science* 2016, 71, 246–249. <https://doi.org/10.1016/j.jcs.2016.09.006>.
- Hell, J.; Prückler, M.; Danner, L.; Henniges, U.; Apprich, S.; Rosenau, T.; Kneifel, W.; Böhmendorfer, S. A Comparison between Near-Infrared (NIR) and Mid-Infrared (ATR-FTIR) Spectroscopy for the Multivariate Determination of Compositional Properties in Wheat Bran Samples. *Food Control* 2016, 60, 365–369. <https://doi.org/10.1016/j.foodcont.2015.08.003>.
- Kim, Y.; Himmelsbach, D. S.; Kays, S. E. ATR-Fourier Transform Mid-Infrared Spectroscopy for Determination of Trans Fatty Acids in Ground Cereal Products without Oil Extraction. *J. Agric. Food Chem.* 2007, 55 (11), 4327–4333. <https://doi.org/10.1021/jf063729l>.
- Lohumi, S.; Lee, S.; Lee, H.; Cho, B.-K. A Review of Vibrational Spectroscopic Techniques for the Detection of Food Authenticity and Adulteration. *Trends in Food Science & Technology* 2015, 46 (1), 85–98. <https://doi.org/10.1016/j.tifs.2015.08.003>.
- Rodríguez, S. D.; Rolandelli, G.; Buera, M. P. Detection of Quinoa Flour Adulteration by Means of FT-MIR Spectroscopy Combined with Chemometric Methods. *Food Chemistry* 2019, 274, 392–401. <https://doi.org/10.1016/j.foodchem.2018.08.140>.
- Shotts, M.-L.; Plans Pujolras, M.; Rossell, C.; Rodríguez-Saona, L. Authentication of Indigenous Flours (Quinoa, Amaranth and Kañiwa) from the Andean Region Using a Portable ATR-Infrared Device in Combination with Pattern Recognition Analysis. *Journal of Cereal Science* 2018, 82, 65–72. <https://doi.org/10.1016/j.jcs.2018.04.005>.

Suchowilska, E.; Kandler, W.; Wiwart, M.; Krska, R. Fourier Transform Infrared - Attenuated Total Reflection for Wheat Grain. *International Agrophysics* 2012, 26 (2), 207–210. <https://doi.org/10.2478/v10247-012-0030-x>.

Sujka, K.; Koczoń, P.; Ceglińska, A.; Reder, M.; Cierniewska-Żytkiewicz, H. The Application of FT-IR Spectroscopy for Quality Control of Flours Obtained from Polish Producers <https://www.hindawi.com/journals/jamc/2017/4315678/> (accessed Nov 4, 2020). <https://doi.org/10.1155/2017/4315678>.

Suchowilska et al. - 2012 - Fourier Transform Infrared - Attenuated Total Refl.Pdf.

Rapid Determination of Carbohydrates, Ash, and Extractives Contents of Straw Using Attenuated Total Reflectance Fourier Transform Mid-Infrared Spectroscopy | *Journal of Agricultural and Food Chemistry* <https://pubs.acs.org/doi/10.1021/jf200078h> (accessed Oct 21, 2020).

3 PRINCIPLE

Mid-infrared spectroscopy measures the absorption of radiation in the frequency range from 4000 to 400 cm^{-1} . The absorption involves transitions between vibrational energy states and rotational substates of the molecule. A selection rule applies to these transitions, and absorption of infrared light can occur only if the vibration causes a change in the dipole moment of molecules. It is possible to assign absorptions to specific functional groups, making IR spectroscopy very useful in structural elucidation. Because the intensity of the absorption is proportional to the concentration of the absorbing species, quantitative analysis is also possible (Kačuráková and Wilson, 2001).

MIRS analysis involves some steps. First, the spectrum of the air is collected as a reference spectrum (background) and then MIR spectra are collected and corrected with reference spectrum. The pre-treatment and multivariate analysis can be performed by using any chemometrics software.

4 PROCEDURE

4.1 MIRS measurements

Eighty-nine samples of cell walls from 37 different varieties (2 or 3 roots per genotype) were provided by CIAT (Colombia) in September 2019. The protocols of preparation of flours and cell walls are detailed in the WP2 SOPs “WP2_SOP _flours starch_finallyzed” and “WP2_SOP_CellWalls_tt” respectively.

Spectra are collected using Thermo Scientific™ Nicolet™ iS50R Research FTIR Spectrometers in the 4000-400 cm^{-1} region and resolution of 0.09 cm^{-1} (fig.1). The acquisition system is controlled by Thermo Scientific™ OMNIC™ Series Software. This software provide not only the acquisition and saving of spectra but also to set the parameters of acquisition such as spectral range, number of scans, offset, gain, integration time, resolution, spectral path.



Figure 1 : Mid infrared spectrometer system

Background spectrum of the air is collected before spectra acquisition. This spectra is used as reference for all measured spectra. Therefore, by using a small spatula, milled cell walls samples (~0.1 g) are placed directly on the surface of the diamond ATR crystal by pressing the powder sample onto the crystal using a pressure clamp and spectrum was collected (fig.2). Three Spectra were collected per sample.



Figure 2 : MIRS measurement of cassava cell walls

4.2 Repeatability

Eight spectra of the same flour sample were acquired in order to determine the optimum number of MIRS measurements (replicates) required to be representative of the sample. Then, the calculation of the root mean square error (RMS) is performed for the spectral range: 1400-800 cm^{-1} (fig.3). The individuals RMS are first calculated for the 8 measurements and then all the RMS are calculated for 5 combinations of 2 measurements randomly selected. The objective is to verify if the RMS value of any random of two samples put together remains of the same order than the RMS observed for 8 replicates.




Figure 3: eight spectra measured on the same sample of cassava flour

The average RMS value obtained for the 8 replicates is 1607 uabs, ranged from 277 to 2493 μ abs with a standard deviation of 1101 μ abs. All of these values are within the confidence interval of $1607 \pm 3SD$.

The RMS obtained for 5 combinations of two replicates (2-3, 1-2, 3-5, 3-7 and 6-8) are 457, 1066, 710, 1352 and 614 μ abs respectively. The RMS values between 2 replicates, randomly selected, were lower than the RMS mean for all the replicates. These results indicate that two replicates of the sample result to a spectral dispersion (variability) similar to the one obtained with 8 repetitions.

5 CRITICAL POINTS OR NOTE ON THE PROCEDURE

- To avoid the moisture absorption, store flours or cell walls samples in a dry place and closed plastic vials.
- To avoid the light diffusion on the spectra, place the sample on all the surface of the diamond ATR crystal.
- Clean the crystal and pressure clamp by using ethanol and paper towel after each measurement.

6 APPENDICES

6.1 Appendix 1: Procedure of extraction of cell walls from cassava flour



Weight the cassava flour



Add 60 mL water



Homogenize



Add alpha-amylase



Add 10 mL Milli-Q water




Clean Turrax's dispersion rod



Incubate 70°C, 17 hrs



Adjust pH



Add amyloglucosidase



Pour sample



Centrifuge

Iodine test



Take out supernatant

Incubate 60°C, 3 hrs



Add 10 mL of Ethanol



Take out supernatant



Centrifuge



Vortex

REPEAT PREVIOUS 4
STEPS WITH EACH
ETHANOL SOLUTION



Add liquid nitrogen



Freeze dry



Grind



Institute: Cirad – UMR QualiSud
Address: C/O Cathy Méjean, TA-B95/15 - 73 rue Jean-François Breton - 34398 Montpellier Cedex 5 - France
Tel: +33 4 67 61 44 31
Email: rtbfoodspmu@cirad.fr
Website: <https://rtbfoods.cirad.fr/>